

Antibodies Product Datasheet

Product Name:	CTG1B Antibody
Catalog Number:	C43327
Reactivity:	Human, Rat, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human CTG1B
Concentration:	1 mg/ml
Gene Name:	CTAG1A CTAG CTAG1 ESO1 LAGE2 LAGE2A; CTAG1B LAGE2B
Protein Name:	Cancer/testis antigen 1 (Autoimmunogenic cancer/testis antigen NY-ESO-1) (Cancer/testis antigen 6.1) (CT6.1) (L antigen family member 2) (LAGE-2)
Other Name:	
Human Swiss Prot No.:	P78358
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human CTG1B
Storage and Stability:	-20°C/1 year