

Antibodies Product Datasheet

Product Name:	CTGE1 Antibody
Catalog Number:	C43328
Reactivity:	Human, Rat, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human CTGE1
Concentration:	1 mg/ml
Gene Name:	CTAGE1 CTAGE2
Protein Name:	Cutaneous T-cell lymphoma-associated antigen 1 (Protein cTAGE-1) (Cancer/testis antigen 21.1) (CT21.1)
Other Name:	
Human Swiss Prot No.:	Q9HC47
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human CTGE1
Storage and Stability:	-20°C/1 year