

Antibodies Product Datasheet

Product Name:	CX6A2 Antibody
Catalog Number:	C43352
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 10-90
Concentration:	1 mg/ml
Gene Name:	COX6A2 COX6A COX6AH
Protein Name:	Cytochrome c oxidase subunit 6A2, mitochondrial (Cytochrome c oxidase polypeptide VIa-heart) (COXVIAH) (Cytochrome c oxidase subunit VIA-muscle) (COX VIa-M)
Other Name:	
Human Swiss Prot No.:	Q02221
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	CX6A2 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year