

Antibodies Product Datasheet

Product Name:	CYH1 Antibody
Catalog Number:	C43361
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from part region of human protein
Concentration:	1 mg/ml
Gene Name:	CYTH1 D17S811E PSCD1
Protein Name:	Cytohesin-1 (PH, SEC7 and coiled-coil domain-containing protein 1) (SEC7 homolog B2-1)
Other Name:	
Human Swiss Prot No.:	Q15438
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	CYH1 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year