

Antibodies Product Datasheet

Product Name:	CYTSB Antibody
Catalog Number:	C43371
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human CYTSB
Concentration:	1 mg/ml
Gene Name:	SPECC1 CYTSB NSP5
Protein Name:	Cytospin-B (Nuclear structure protein 5) (NSP5) (Sperm antigen HCMOGT-1) (Sperm antigen with calponin homology and coiled-coil domains 1)
Other Name:	
Human Swiss Prot No.:	Q5M775
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human CYTSB
Storage and Stability:	-20°C/1 year