

Antibodies Product Datasheet

Product Name:	D107A Antibody
Catalog Number:	C43373
Reactivity:	Human, Rat, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human D107A
Concentration:	1 mg/ml
Gene Name:	DEFB107A DEFB107 DEFB7; DEFB107B
Protein Name:	Beta-defensin 107 (Beta-defensin 7) (BD-7) (DEFB-7) (Defensin, beta 107)
Other Name:	
Human Swiss Prot No.:	Q8IZN7
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human D107A
Storage and Stability:	-20°C/1 year