

Antibodies Product Datasheet

Product Name:	D108A Antibody
Catalog Number:	C43374
Reactivity:	Human, Rat, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human D108A
Concentration:	1 mg/ml
Gene Name:	DEFB108P1 DEFB108A; DEFB108P2
Protein Name:	Putative beta-defensin 108A (Defensin, beta 108A) (Putative beta-defensin 108B pseudogene 1/2)
Other Name:	
Human Swiss Prot No.:	A8MXU0
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human D108A
Storage and Stability:	-20°C/1 year