

Antibodies Product Datasheet

Product Name:	DACT2 Antibody
Catalog Number:	C43376
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human DACT2
Concentration:	1 mg/ml
Gene Name:	DACT2 C6orf116 PP13671
Protein Name:	DACT2
Other Name:	Dapper homolog 2 (Dapper antagonist of catenin 2)
Human Swiss Prot No.:	Q5SW24
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.355% sodium azide.
Specificity:	This antibody detects endogenous levels of DACT2 at Human/Mouse
Storage and Stability:	-20°C/1 year