

Antibodies Product Datasheet

Product Name:	DC8L2 Antibody
Catalog Number:	C43390
Reactivity:	Human, Rat, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human DC8L2
Concentration:	1 mg/ml
Gene Name:	DCAF8L2 WDR42C
Protein Name:	DDB1- and CUL4-associated factor 8-like protein 2 (WD repeat-containing protein 42C)
Other Name:	
Human Swiss Prot No.:	P0C7V8
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human DC8L2
Storage and Stability:	-20°C/1 year