

Antibodies Product Datasheet

Product Name:	DCBD2 Antibody
Catalog Number:	C43395
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human DCBD2
Concentration:	1 mg/ml
Gene Name:	DCBLD2 CLCP1 ESDN
Protein Name:	DCBD2
Other Name:	Discoidin, CUB and LCCL domain-containing protein 2 (CUB, LCCL and coagulation factor V/VIII-homology domains protein 1) (Endothelial and smooth muscle cell-derived neuropilin-like protein)
Human Swiss Prot No.:	Q96PD2
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.359% sodium azide.
Specificity:	This antibody detects endogenous levels of DCBD2 at Human/Mouse/Rat
Storage and Stability:	-20°C/1 year