

Antibodies Product Datasheet

Product Name:	DCST2 Antibody
Catalog Number:	C43400
Reactivity:	Human, Rat, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human DCST2
Concentration:	1 mg/ml
Gene Name:	DCST2
Protein Name:	DCST2
Other Name:	DC-STAMP domain-containing protein 2
Human Swiss Prot No.:	Q5T1A1
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.347% sodium azide.
Specificity:	This antibody detects endogenous levels of DCST2 at Human
Storage and Stability:	-20°C/1 year