

Antibodies Product Datasheet

Product Name:	DDAH2 Antibody
Catalog Number:	C43403
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human DDAH2
Concentration:	1 mg/ml
Gene Name:	DDAH2 DDAH G6A NG30
Protein Name:	N(G),N(G)-dimethylarginine dimethylaminohydrolase 2 (DDAH-2) (Dimethylarginine dimethylaminohydrolase 2) (EC 3.5.3.18) (DDAHII) (Dimethylargininase-2) (Protein G6a) (S-phase protein)
Other Name:	
Human Swiss Prot No.:	O95865
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human DDAH2
Storage and Stability:	-20°C/1 year