

Antibodies Product Datasheet

Product Name:	DDX11 Antibody
Catalog Number:	C43407
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human DDX11
Concentration:	1 mg/ml
Gene Name:	DDX11 CHL1 CHLR1 KRG2
Protein Name:	DDX11
Other Name:	Probable ATP-dependent RNA helicase DDX11 (EC 3.6.4.13) (CHL1-related protein 1) (hCHLR1) (DEAD/H box protein 11) (Keratinocyte growth factor-regulated gene 2 protein) (KRG-2)
Human Swiss Prot No.:	Q96FC9
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.281% sodium azide.
Specificity:	This antibody detects endogenous levels of DDX11 at Human/Mouse
Storage and Stability:	-20°C/1 year