

## Antibodies Product Datasheet

|                               |                                                                                                                   |
|-------------------------------|-------------------------------------------------------------------------------------------------------------------|
| <b>Product Name:</b>          | DDX12 Antibody                                                                                                    |
| <b>Catalog Number:</b>        | C43408                                                                                                            |
| <b>Reactivity:</b>            | Human, Rat, Mouse                                                                                                 |
| <b>Applications:</b>          | WB                                                                                                                |
| <b>Dilutions:</b>             | WB 1:500-2000                                                                                                     |
| <b>Modification:</b>          | Unmodified/Total                                                                                                  |
| <b>Source:</b>                | Polyclonal Rabbit                                                                                                 |
| <b>Purification:</b>          | The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.         |
| <b>Immunogen:</b>             | Synthesized peptide derived from human DDX12                                                                      |
| <b>Concentration:</b>         | 1 mg/ml                                                                                                           |
| <b>Gene Name:</b>             | DDX12P CHLR2 DDX12                                                                                                |
| <b>Protein Name:</b>          | DDX12                                                                                                             |
| <b>Other Name:</b>            | Putative ATP-dependent RNA helicase DDX12 (EC 3.6.4.13) (CHL1-related protein 2) (hCHLR2) (DEAD/H box protein 12) |
| <b>Human Swiss Prot No.:</b>  | Q92771                                                                                                            |
| <b>Formulation:</b>           | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.208% sodium azide.                                          |
| <b>Specificity:</b>           | This antibody detects endogenous levels of DDX12 at Human                                                         |
| <b>Storage and Stability:</b> | -20°C/1 year                                                                                                      |