

Antibodies Product Datasheet

Product Name:	DDX42 Antibody
Catalog Number:	C43410
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 850-930
Concentration:	1 mg/ml
Gene Name:	DDX42
Protein Name:	ATP-dependent RNA helicase DDX42 (EC 3.6.4.13) (DEAD box protein 42) (RNA helicase-like protein) (RHELP) (RNA helicase-related protein) (RNAHP) (SF3b DEAD box protein) (Splicing factor 3B-associated 125 kDa protein) (SF3b125)
Other Name:	
Human Swiss Prot No.:	Q86XP3
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	DDX42 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year