

Antibodies Product Datasheet

Product Name:	DEAF1 Antibody
Catalog Number:	C43416
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human DEAF1
Concentration:	1 mg/ml
Gene Name:	DEAF1 SPN ZMYND5
Protein Name:	DEAF1
Other Name:	Deformed epidermal autoregulatory factor 1 homolog (Nuclear DEAF-1-related transcriptional regulator) (NUDR) (Suppressin) (Zinc finger MYND domain-containing protein 5)
Human Swiss Prot No.:	O75398
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.125% sodium azide.
Specificity:	This antibody detects endogenous levels of DEAF1 at Human/Mouse/Rat
Storage and Stability:	-20°C/1 year