

Antibodies Product Datasheet

Product Name:	DERL2 Antibody
Catalog Number:	C43424
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from part region of human protein
Concentration:	1 mg/ml
Gene Name:	DERL2 DER2 FLANA CGI-101 SBB153
Protein Name:	Derlin-2 (Degradation in endoplasmic reticulum protein 2) (DERtrin-2) (Der1-like protein 2) (F-LAN-1) (F-LANa)
Other Name:	
Human Swiss Prot No.:	Q9GZP9
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	DERL2 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year