

Antibodies Product Datasheet

Product Name:	DGAT2 Antibody
Catalog Number:	C43430
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from part region of human protein
Concentration:	1 mg/ml
Gene Name:	DGAT2 HMFN1045 UNQ738/PRO1433
Protein Name:	Diacylglycerol O-acyltransferase 2 (EC 2.3.1.20) (Acyl-CoA retinol O-fatty-acyltransferase) (ARAT) (Retinol O-fatty-acyltransferase) (EC 2.3.1.76) (Diglyceride acyltransferase 2)
Other Name:	
Human Swiss Prot No.:	Q96PD7
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	DGAT2 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year