

Antibodies Product Datasheet

Product Name:	DHB12 Antibody
Catalog Number:	C43434
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from part region of human protein
Concentration:	1 mg/ml
Gene Name:	HSD17B12
Protein Name:	Estradiol 17-beta-dehydrogenase 12 (EC 1.1.1.62) (17-beta-hydroxysteroid dehydrogenase 12) (17-beta-HSD 12) (3-ketoacyl-CoA reductase) (KAR) (EC 1.3.1.-)
Other Name:	
Human Swiss Prot No.:	Q53GQ0
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	DHB12 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year