

Antibodies Product Datasheet

Product Name:	DHB14 Antibody
Catalog Number:	C43435
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from part region of human protein
Concentration:	1 mg/ml
Gene Name:	HSD17B14 DHRS10 SDR3 UNQ502/PRO474
Protein Name:	17-beta-hydroxysteroid dehydrogenase 14 (17-beta-HSD 14) (EC 1.1.1.-) (17-beta-hydroxysteroid dehydrogenase DHRS10) (Dehydrogenase/reductase SDR family member 10) (Retinal short-chain dehydrogenase/reductase retSDR3)
Other Name:	
Human Swiss Prot No.:	Q9BPX1
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	DHB14 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year