

Antibodies Product Datasheet

Product Name:	DHB2 Antibody
Catalog Number:	C43436
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from part region of human protein
Concentration:	1 mg/ml
Gene Name:	HSD17B2 EDH17B2
Protein Name:	Estradiol 17-beta-dehydrogenase 2 (EC 1.1.1.62) (17-beta-hydroxysteroid dehydrogenase type 2) (17-beta-HSD 2) (20 alpha-hydroxysteroid dehydrogenase) (20-alpha-HSD) (E2DH) (Microsomal 17-beta-hydroxysteroid dehydrogenase) (Testosterone 17-beta-dehydrogenase)
Other Name:	
Human Swiss Prot No.:	P37059
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	DHB2 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year