

Antibodies Product Datasheet

Product Name:	DHI2 Antibody
Catalog Number:	C43438
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from part region of human protein
Concentration:	1 mg/ml
Gene Name:	HSD11B2 HSD11K
Protein Name:	Corticosteroid 11-beta-dehydrogenase isozyme 2 (EC 1.1.1.-) (11-beta-hydroxysteroid dehydrogenase type 2) (11-DH2) (11-beta-HSD2) (11-beta-hydroxysteroid dehydrogenase type II) (-HSD11 type II) (NAD-dependent 11-beta-hydroxysteroid dehydrogenase) (11-beta
Other Name:	
Human Swiss Prot No.:	P80365
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	DHI2 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year