

Antibodies Product Datasheet

Product Name:	DHRS9 Antibody
Catalog Number:	C43441
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from part region of human protein
Concentration:	1 mg/ml
Gene Name:	DHRS9 UNQ835/PRO1773
Protein Name:	Dehydrogenase/reductase SDR family member 9 (EC 1.1.-.-) (3-alpha hydroxysteroid dehydrogenase) (3-alpha-HSD) (NADP-dependent retinol dehydrogenase/reductase) (RDH-E2) (RDHL) (Short-chain dehydrogenase/reductase retSDR8) (Tracheobronchial epithelial cell-
Other Name:	
Human Swiss Prot No.:	Q9BPW9
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	DHRS9 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year