

Antibodies Product Datasheet

Product Name:	DHSO Antibody
Catalog Number:	C43442
Reactivity:	Human, Mouse, Rat
Applications:	WB, IP
Dilutions:	WB 1: 500-2000 IP 1:200
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human DHSO
Concentration:	1 mg/ml
Gene Name:	SORD
Protein Name:	DHSO
Other Name:	
Human Swiss Prot No.:	Q00796
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of DHSO at Human/Mouse/Rat
Storage and Stability:	-20°C/1 year