

Antibodies Product Datasheet

Product Name:	DIXC1 Antibody
Catalog Number:	C43458
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human DIXC1
Concentration:	1 mg/ml
Gene Name:	DIXDC1 CCD1 KIAA1735
Protein Name:	Dixin (Coiled-coil protein DIX1) (Coiled-coil-DIX1) (DIX domain-containing protein 1)
Other Name:	
Human Swiss Prot No.:	Q155Q3
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human DIXC1
Storage and Stability:	-20°C/1 year