

Antibodies Product Datasheet

Product Name:	DJB14 Antibody
Catalog Number:	C43459
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human DJB14
Concentration:	1 mg/ml
Gene Name:	DNAJB14 UNQ9427/PRO34683
Protein Name:	DnaJ homolog subfamily B member 14
Other Name:	
Human Swiss Prot No.:	Q8TBM8
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human DJB14
Storage and Stability:	-20°C/1 year