

Antibodies Product Datasheet

Product Name:	DKK4 Antibody
Catalog Number:	C43464
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 150-230
Concentration:	1 mg/ml
Gene Name:	DKK4
Protein Name:	Dickkopf-related protein 4 (Dickkopf-4) (Dkk-4) (hDkk-4) [Cleaved into: Dickkopf-related protein 4 short form]
Other Name:	
Human Swiss Prot No.:	Q9UBT3
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	DKK4 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year