

Antibodies Product Datasheet

Product Name:	DLRB1 Antibody
Catalog Number:	C43472
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from part region of human protein
Concentration:	1 mg/ml
Gene Name:	DYNLRB1 BITH DNCL2A DNLC2A ROBLD1 HSPC162
Protein Name:	Dynein light chain roadblock-type 1 (Bithoraxoid-like protein) (BLP) (Dynein light chain 2A, cytoplasmic) (Dynein-associated protein Km23) (Roadblock domain-containing protein 1)
Other Name:	
Human Swiss Prot No.:	Q9NP97
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	DLRB1 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year