

Antibodies Product Datasheet

Product Name:	DMBT1 Antibody
Catalog Number:	C43475
Reactivity:	Human, Mouse, Rat
Applications:	IHC-p, IF(paraffin section)
Dilutions:	IHC-p 1:50-300
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from part region of human protein
Concentration:	1 mg/ml
Gene Name:	DMBT1 GP340
Protein Name:	Deleted in malignant brain tumors 1 protein (Glycoprotein 340) (Gp-340) (Hensin) (Salivary agglutinin) (SAG) (Surfactant pulmonary-associated D-binding protein)
Other Name:	
Human Swiss Prot No.:	Q9UGM3
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	DMBT1 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year