

Antibodies Product Datasheet

Product Name:	DMRT2 Antibody
Catalog Number:	C43480
Reactivity:	Human, Rat, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human DMRT2
Concentration:	1 mg/ml
Gene Name:	DMRT2 DSXL2
Protein Name:	DMRT2
Other Name:	Doublesex- and mab-3-related transcription factor 2 (Doublesex-like 2 protein) (DSXL-2)
Human Swiss Prot No.:	Q9Y5R5
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.85% sodium azide.
Specificity:	This antibody detects endogenous levels of DMRT2 at Human
Storage and Stability:	-20°C/1 year