

Antibodies Product Datasheet

Product Name:	DNA2 Antibody
Catalog Number:	C43484
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from part region of human protein
Concentration:	1 mg/ml
Gene Name:	DNA2 DNA2L KIAA0083
Protein Name:	DNA replication ATP-dependent helicase/nuclease DNA2 (hDNA2) (DNA replication ATP-dependent helicase-like homolog) [Includes: DNA replication nuclease DNA2 (EC 3.1.-.-); DNA replication ATP-dependent helicase DNA2 (EC 3.6.4.12)]
Other Name:	
Human Swiss Prot No.:	P51530
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	DNA2 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year