

Antibodies Product Datasheet

Product Name:	DND1 Antibody
Catalog Number:	C43485
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human DND1
Concentration:	1 mg/ml
Gene Name:	DND1 RBMS4
Protein Name:	DND1
Other Name:	Dead end protein homolog 1 (RNA-binding motif, single-stranded-interacting protein 4)
Human Swiss Prot No.:	Q8IYX4
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.33% sodium azide.
Specificity:	This antibody detects endogenous levels of DND1 at Human/Mouse
Storage and Stability:	-20°C/1 year