

Antibodies Product Datasheet

Product Name:	DNJB2 Antibody
Catalog Number:	C43488
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 260-340
Concentration:	1 mg/ml
Gene Name:	DNAJB2 HSJ1 HSPF3
Protein Name:	DnaJ homolog subfamily B member 2 (DnaJ protein homolog 1) (Heat shock 40 kDa protein 3) (Heat shock protein J1) (HSJ-1)
Other Name:	
Human Swiss Prot No.:	P25686
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	DNJB2 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year