

Antibodies Product Datasheet

Product Name:	DNJB5 Antibody
Catalog Number:	C43489
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 140-220
Concentration:	1 mg/ml
Gene Name:	DNAJB5 HSC40
Protein Name:	DnaJ homolog subfamily B member 5 (Heat shock protein Hsp40-2) (Heat shock protein Hsp40-3) (Heat shock protein cognate 40) (Hsc40)
Other Name:	
Human Swiss Prot No.:	O75953
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	DNJB5 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year