

Antibodies Product Datasheet

Product Name:	DNJC3 Antibody
Catalog Number:	C43492
Reactivity:	Human, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 80-160
Concentration:	1 mg/ml
Gene Name:	DNAJC3 P58IPK PRKRI
Protein Name:	DnaJ homolog subfamily C member 3 (Endoplasmic reticulum DnaJ protein 6) (ERdj6) (Interferon-induced, double-stranded RNA-activated protein kinase inhibitor) (Protein kinase inhibitor of 58 kDa) (Protein kinase inhibitor p58)
Other Name:	
Human Swiss Prot No.:	Q13217
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	DNJC3 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year