

Antibodies Product Datasheet

Product Name:	DNJC8 Antibody
Catalog Number:	C43493
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human DNJC8
Concentration:	1 mg/ml
Gene Name:	DNAJC8 SPF31 HSPC315 HSPC331
Protein Name:	DnaJ homolog subfamily C member 8 (Splicing protein spf31)
Other Name:	
Human Swiss Prot No.:	O75937
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human DNJC8
Storage and Stability:	-20°C/1 year