

Antibodies Product Datasheet

Product Name:	DPH3B Antibody
Catalog Number:	C43504
Reactivity:	Human, Rat, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human DPH3B
Concentration:	1 mg/ml
Gene Name:	DPH3P1 C20orf143 DPH3B ZCSL1
Protein Name:	Putative DPH3 homolog B (CSL-type zinc finger-containing protein 1)
Other Name:	
Human Swiss Prot No.:	Q9H4G8
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human DPH3B
Storage and Stability:	-20°C/1 year