

Antibodies Product Datasheet

Product Name:	DPP3 Antibody
Catalog Number:	C43510
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human DPP3
Concentration:	1 mg/ml
Gene Name:	DPP3
Protein Name:	Dipeptidyl peptidase 3 (EC 3.4.14.4) (Dipeptidyl aminopeptidase III) (Dipeptidyl arylamidase III) (Dipeptidyl peptidase III) (DPP III)
Other Name:	
Human Swiss Prot No.:	Q9NY33
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human DPP3
Storage and Stability:	-20°C/1 year