

Antibodies Product Datasheet

Product Name:	DPYL1 Antibody
Catalog Number:	C43514
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 440-520
Concentration:	1 mg/ml
Gene Name:	CRMP1 DPYSL1 ULIP3
Protein Name:	Dihydropyrimidinase-related protein 1 (DRP-1) (Collapsin response mediator protein 1) (CRMP-1) (Unc-33-like phosphoprotein 3) (ULIP-3)
Other Name:	
Human Swiss Prot No.:	Q14194
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	DPYL1 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year