

Antibodies Product Datasheet

Product Name:	DPYL3 Antibody
Catalog Number:	C43515
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human DPYL3
Concentration:	1 mg/ml
Gene Name:	DPYSL3 CRMP4 DRP3 ULIP ULIP1
Protein Name:	DPYL3
Other Name:	Dihydropyrimidinase-related protein 3 (DRP-3) (Collapsin response mediator protein 4) (CRMP-4) (Unc-33-like phosphoprotein 1) (ULIP-1)
Human Swiss Prot No.:	Q14195
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.179% sodium azide.
Specificity:	This antibody detects endogenous levels of DPYL3 at Human/Mouse/Rat
Storage and Stability:	-20°C/1 year