

Antibodies Product Datasheet

Product Name:	DPYL5 Antibody
Catalog Number:	C43516
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human DPYL5
Concentration:	1 mg/ml
Gene Name:	DPYSL5 CRMP5 ULIP6
Protein Name:	DPYL5
Other Name:	Dihydropyrimidinase-related protein 5 (DRP-5) (CRMP3-associated molecule) (CRAM) (Collapsin response mediator protein 5) (CRMP-5) (UNC33-like phosphoprotein 6) (ULIP-6)
Human Swiss Prot No.:	Q9BPU6
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.117% sodium azide.
Specificity:	This antibody detects endogenous levels of DPYL5 at Human/Mouse/Rat
Storage and Stability:	-20°C/1 year