

Antibodies Product Datasheet

Product Name:	DSC2 Antibody
Catalog Number:	C43526
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 150-230
Concentration:	1 mg/ml
Gene Name:	DSC2 CDHF2 DSC3
Protein Name:	Desmocollin-2 (Cadherin family member 2) (Desmocollin-3) (Desmosomal glycoprotein II) (Desmosomal glycoprotein III)
Other Name:	
Human Swiss Prot No.:	Q02487
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	DSC2 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year