

Antibodies Product Datasheet

Product Name:	DSE Antibody
Catalog Number:	C43532
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human DSE
Concentration:	1 mg/ml
Gene Name:	DSE SART2
Protein Name:	DSE
Other Name:	Dermatan-sulfate epimerase (DS epimerase) (EC 5.1.3.19) (Chondroitin-glucuronate 5-epimerase) (Squamous cell carcinoma antigen recognized by T-cells 2) (SART-2)
Human Swiss Prot No.:	Q9UL01
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.249% sodium azide.
Specificity:	This antibody detects endogenous levels of DSE at Human/Mouse
Storage and Stability:	-20°C/1 year