

Antibodies Product Datasheet

Product Name:	DTL Antibody
Catalog Number:	C43535
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from part region of human protein
Concentration:	1 mg/ml
Gene Name:	DTL CDT2 CDW1 DCAF2 L2DTL RAMP
Protein Name:	Denticleless protein homolog (DDB1- and CUL4-associated factor 2) (Lethal(2) denticleless protein homolog) (Retinoic acid-regulated nuclear matrix-associated protein)
Other Name:	
Human Swiss Prot No.:	Q9NZJ0
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	DTL Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year