

Antibodies Product Datasheet

Product Name:	DUS14 Antibody
Catalog Number:	C43538
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 120-200
Concentration:	1 mg/ml
Gene Name:	DUSP14 MKP6
Protein Name:	Dual specificity protein phosphatase 14 (EC 3.1.3.16) (EC 3.1.3.48) (MKP-1-like protein tyrosine phosphatase) (MKP-L) (Mitogen-activated protein kinase phosphatase 6) (MAP kinase phosphatase 6) (MKP-6)
Other Name:	
Human Swiss Prot No.:	O95147
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	DUS14 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year