

Antibodies Product Datasheet

Product Name:	DUS23 Antibody
Catalog Number:	C43541
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human DUS23
Concentration:	1 mg/ml
Gene Name:	DUSP23 LDP3 VH2
Protein Name:	Dual specificity protein phosphatase 23 (EC 3.1.3.16) (EC 3.1.3.48) (Low molecular mass dual specificity phosphatase 3) (LDP-3) (VH1-like phosphatase Z)
Other Name:	
Human Swiss Prot No.:	Q9BVJ7
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human DUS23
Storage and Stability:	-20°C/1 year