

Antibodies Product Datasheet

Product Name:	DYH5 Antibody
Catalog Number:	C43548
Reactivity:	Human, Mouse
Applications:	IHC-p, IF(paraffin section)
Dilutions:	IHC-p 1: 50-200
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human DYH5
Concentration:	1 mg/ml
Gene Name:	DNAH5 DNAHC5 HL1 KIAA1603
Protein Name:	DYH5
Other Name:	
Human Swiss Prot No.:	Q8TE73
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of DYH5 at Human/Mouse
Storage and Stability:	-20°C/1 year