

Antibodies Product Datasheet

Product Name:	DYHC1 Antibody
Catalog Number:	C43550
Reactivity:	Human, Mouse, Rat
Applications:	IHC-p, IF(paraffin section)
Dilutions:	IHC-p 1:50-300
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from part region of human protein
Concentration:	1 mg/ml
Gene Name:	DYNC1H1 DHC1 DNCH1 DNCL DNECL DYHC KIAA0325
Protein Name:	Cytoplasmic dynein 1 heavy chain 1 (Cytoplasmic dynein heavy chain 1) (Dynein heavy chain, cytosolic)
Other Name:	
Human Swiss Prot No.:	Q14204
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	DYHC1 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year