

Antibodies Product Datasheet

Product Name:	DZAN1 Antibody
Catalog Number:	C43559
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human DZAN1
Concentration:	1 mg/ml
Gene Name:	DZANK1 C20orf12 C20orf84
Protein Name:	Double zinc ribbon and ankyrin repeat-containing protein 1
Other Name:	
Human Swiss Prot No.:	Q9NVP4
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human DZAN1
Storage and Stability:	-20°C/1 year