

Antibodies Product Datasheet

Product Name:	DZI1L Antibody
Catalog Number:	C43560
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human DZI1L
Concentration:	1 mg/ml
Gene Name:	DZIP1L
Protein Name:	DZI1L
Other Name:	Zinc finger protein DZIP1L (DAZ-interacting protein 1-like protein)
Human Swiss Prot No.:	Q8IYY4
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.319% sodium azide.
Specificity:	This antibody detects endogenous levels of DZI1L at Human/Mouse/Rat
Storage and Stability:	-20°C/1 year